

Item No.	Description	Size	Packaging	
			Inner	Outer
Safety Blood Collection/Infusion Sets (sterile)				
450090	Safety Blood Collection/Infusion Set	tubing length: 19cm 25G x 3/4"	50 pcs.	1000 pcs.
450091	Safety Blood Collection/Infusion Set	tubing length: 19cm 21G x 3/4"	50 pcs.	1000 pcs.
450092	Safety Blood Collection/Infusion Set	tubing length: 19cm 23G x 3/4"	50 pcs.	1000 pcs.
Safety Blood Collection Sets + Luer Adapters (sterile)				
450083	Safety Blood Collection Set + Luer Adapter	tubing length: 10cm 21G x 3/4"	50 pcs.	1000 pcs.
450081	Safety Blood Collection Set + Luer Adapter	tubing length: 19cm 21G x 3/4"	50 pcs.	1000 pcs.
450095	Safety Blood Collection Set + Luer Adapter	tubing length: 30cm 21G x 3/4"	50 pcs.	1000 pcs.
450084	Safety Blood Collection Set + Luer Adapter	tubing length: 10cm 23G x 3/4"	50 pcs.	1000 pcs.
450082	Safety Blood Collection Set + Luer Adapter	tubing length: 19cm 23G x 3/4"	50 pcs.	1000 pcs.
450096	Safety Blood Collection Set + Luer Adapter	tubing length: 30cm 23G x 3/4"	50 pcs.	1000 pcs.
450099	Safety Blood Collection Set + Luer Adapter	tubing length: 30cm 25G x 3/4"	50 pcs.	1000 pcs.
Safety Blood Collection Sets + Holders (sterile)				
450087	Safety Blood Collection Set + Holder	tubing length: 10cm 21G x 3/4"	24 pcs.	240 pcs.
450085	Safety Blood Collection Set + Holder	tubing length: 19cm 21G x 3/4"	24 pcs.	240 pcs.
450086	Safety Blood Collection Set + Holder	tubing length: 19cm 23G x 3/4"	24 pcs.	240 pcs.
450080	Safety Blood Collection Set + Holder	tubing length: 30cm 25G x 3/4"	24 pcs.	240 pcs.
Safety Blood Collection Sets + Blood Culture Holders (sterile)				
450182	Safety Blood Collection Set + Blood Culture Holder	tubing length: 19cm 21G x 3/4"	24 pcs.	240 pcs.
450183	Safety Blood Collection Set + Blood Culture Holder	tubing length: 19cm 23G x 3/4"	24 pcs.	240 pcs.
450184	Safety Blood Collection Set + Blood Culture Holder	tubing length: 30cm 21G x 3/4"	24 pcs.	240 pcs.
450185	Safety Blood Collection Set + Blood Culture Holder	tubing length: 30cm 23G x 3/4"	24 pcs.	240 pcs.

These devices have no components made of natural dry rubber.

Your Power for Health



Your Power for Health



For further information, please visit our website www.gbo.com/preanalytics or contact us:

Austria (Main Office)
Greiner Bio-One GmbH
Phone (+43) 75 83 67 91-0
Fax (+43) 75 83 63 18
E-Mail office@at.gbo.com

France
Greiner Bio-One SAS
Phone (+33) 1 69 86 25 25
Fax (+33) 1 69 86 25 92
E-Mail office@fr.gbo.com

Hungary
Greiner Bio-One Hungary Kft.
Phone (+36) 96 21 30 88
Fax (+36) 96 21 31 98
E-Mail office@hu.gbo.com

Thailand
Greiner Bio-One Thailand Ltd
Phone (+66) 3 84 65 63-30
Fax (+66) 3 84 65 63-6
E-Mail office@th.gbo.com

Brazil
Greiner Bio-One Brasil
Phone (+55) 19 34 68 96 00
Fax (+55) 19 34 68 96 21
E-Mail office@br.gbo.com

Germany
Greiner Bio-One GmbH
Phone (+49) 2 01 86 18 6-0
Fax (+49) 2 01 86 18 61-2
E-Mail office@de.gbo.com

Spain
VACUETTE Espana S.A.
Phone (+34) 91 652 77 07
Fax (+34) 91 652 33 35
E-Mail info@vacuette.es

The Netherlands
Greiner Bio-One B.V.
Phone (+31) 1 72 42 09 00
Fax (+31) 1 72 44 38 01
E-Mail info@nl.gbo.com

Egypt
Greiner Bio-One Middle East
Phone (+20) 26 21 87 06
Fax (+20) 26 21 87 08
E-Mail hisham.ouda@gbo.com

Great Britain
Greiner Bio-One Ltd.
Phone (+44) 14 53 82 52 55
Fax (+44) 14 53 82 62 66
E-Mail info@uk.gbo.com

Switzerland
Greiner Bio-One VACUETTE
Schweiz GmbH
Phone (+41) 7 12 28 55 22
Fax (+41) 7 12 28 55 21
E-Mail office@ch.gbo.com

USA
Greiner Bio-One North America Inc.
Phone (+1) 70 42 61 78 00
Fax (+1) 70 42 61 78 99
E-Mail office@us.gbo.com



VACUETTE®
Safety Blood
Collection Sets

Our Innovations
for Your Safety

990157 rev. 00, 01.2011 e

VACUETTE® Safety Blood Collection Sets

Specially developed for blood collection from patients with difficult vein conditions. Protect yourself and your employees from the risk of needlestick injury. Every needlestick injury is one too many.



Activation

Complete activation of the safety mechanism is clearly indicated by an audible "click".

Safety Mechanism

Simple activation of the safety mechanism via lightly pressing in both sides of the hub.

Transparent View Window

The transparent view window provides clear vein entry indication.

Especially Safe

Particularly safe due to activation of the safety mechanism while the needle is still in the vein.

More Flexibility


Assorted tube lengths offer more flexibility for blood collection.

Optional

use as an infusion set.

VACUETTE® Safety Blood Collection Sets Combination Products

As a practical combination product, the VACUETTE® Safety Blood Collection Set is also available pre-assembled with:



Item No. 450081-84,
Item No. 450095-96,
Item No. 450099

VACUETTE® Safety Blood Collection Set + Luer Adapter

- Pre-attached set with Luer adapter in sterile blister pack.
- Handling is optimised.



Item No. 450085-87,
Item No. 450080

VACUETTE® Safety Blood Collection Set + Holder

- Pre-attached set with holder in sterile blister pack.
- No threading necessary.
- A practical, ready-to-use solution that saves time.




Item No. 450182-85

VACUETTE® Safety Blood Collection Set + Blood Culture Holder

- Pre-attached set with Blood Culture Holder in sterile blister pack.
- Has been designed for use with most commonly used blood culture bottles and VACUETTE® blood collection tubes in one procedure.
- The set is sterile and pre-assembled for convenience, cost effectiveness and ensures your patient only has to undergo one venipuncture for all procedures.

Instructions for Use *



- Remove protective cap




- 2 Venipuncture



- 3 Insert the tube



- 4 Remove the last tube



- 5 Activate the safety mechanism



- 6 Remove the Safety Blood Collection Set



- 7 Activated Safety Blood Collection Set



- 8 Disposal

* For the complete blood collection procedure please refer to the Instructions for Use.